ORBIT GALAXY® G2 Microcoil Delivery System XTRASOFT® ORBIT GALAXY® G2 Microcoil

Medos International SÀRL
Chemin-Blanc 38, 2400
Le Locle, Switzerland
☎ 41 32 933 83 00

Codman & Shurtleff, Inc.
325 Paramount Drive
Raynham, MA 02767-0350 USA
☎ 800 225 0460
☎ +1 508 828 3000

CE
0086

LCN 55113-001/C

© 2011–2013 Codman & Shurtleff, Inc.

♻️ Revised 06/13

ENGLISH

IMPORTANT INFORMATION

Please Read Before Use

ORBIT GALAXY® G2 Microcoil Delivery System XTRASOFT® ORBIT GALAXY® G2 Microcoil

STERILE R


Rx Only

SYSTEM DESCRIPTION

The ORBIT GALAXY G2 Microcoil Delivery System consists of three components, a Microcoil System, a connecting cable, and a Detachment Control Box (DCB). Each component is sold separately.

As shown in Figure 1, the Microcoil System consists of an embolic microcoil attached to a device positioning unit (DPU) wire. It is covered by an introducer sheath system. The microcoil's shape can be complex or helical. The DPU wire is a variable stiffness pushing system and has a radiopaque marker band located three (3) cm from its distal end. The introducer sheath system has three main components: an introducer tip, a translucent introducer body, and a re-sheathing tool.

The connecting cable delivers the energy needed to detach the coil from the Microcoil System's detachment zone. The connecting cable is connected between the Microcoil System's hub connector on the DPU Wire and the output connector on the DCB. Its length is approximately 5–6 ft (1.5–1.8 m). The connecting cable may be one of two types: one with a remote detach button (the EnPOWER® Control Cable) catalog no. ECB000182-00, or one without a detach button catalog no. CCB000157-00 (the standard connecting cable).

The Detachment Control Box (DCB) provides the energy necessary to allow for a thermo-mechanical detachment of the microcoil from the DPU wire. The blue EnPOWER Detachment Control Box works with the EnPOWER Control Cable or with the standard connecting cable. Both cables are available separately from the microcoil delivery system.

INDICATIONS FOR USE

The XTRASOFT ORBIT GALAXY G2 Microcoil Delivery System is intended for endovascular embolization of intracranial aneurysms.

WARNINGS

- The ORBIT GALAXY G2 Microcoil Delivery System should only be used as a complete system. A complete system requires a Detachment Control Box, a connecting cable, and Microcoil System. DO NOT SUBSTITUTE any components or devices from other manufacturers with the Codman Detachment Control Box or injury to the patient or user could result.
- Explosion Hazard: The Detachment Control Box is not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide or oxygen enriched atmospheres.**

PRECAUTIONS

The Microcoil System shall be used only by physicians trained in interventional neuroradiology and in all aspects of the system.

Do not attempt to sterilize any component of the Microcoil Delivery System.

The Detachment Control Box (DCB) is provided NON-STERILE. Unless placed in a sterile sleeve, it should remain outside of the sterile field during use. The DCB is the only reusable portion of the Microcoil Delivery System. Upon receipt of the DCB, inspect the shipping carton and the unit for signs of damage. If unit presents with physical damage upon shipment, immediately notify Codman & Shurtleff, Inc. for replacement.

Verify the functionality of the Microcoil Delivery System before proceeding with microcoil placement. Refer to the section Microcoil System Preparation on how to verify the functionality.

ADVERSE EVENTS

Possible adverse events include, but are not limited to, the following: hematoma at site of entry, vessel perforation, infection, emboli, hemorrhage, ischemia or vasospasm, neurological deficits including stroke, and possibly death.

HOW SUPPLIED

2

The ORBIT GALAXY G2 Microcoil System and cables are intended for SINGLE USE ONLY; DO NOT RESTERILIZE. Discard the Connecting Cable after one procedure. Structural integrity and/or function may be impaired through reuse or cleaning. Cables are extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused.

Microcoils are permanent implants and are NOT reusable. Discard the microcoil system if it is retrieved from the patient without being implanted. Structural integrity and/or function may be impaired through reuse or cleaning. Microcoils are extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused.

Codman & Shurtleff will not be responsible for any product that is resterilized, nor accept for credit or exchange any product that has been opened but not used.

As long as the inner unit is not opened or damaged, the Microcoil System and Connecting Cable are sterile. The Microcoil System is nonpyrogenic.

REQUIRED ADDITIONAL MATERIALS

The following items are either supplied with the system or sold separately and must be on-hand prior to commencing a procedure.

- Femoral sheath
- Guiding catheter, 5 to 7 Fr
- Infusion microcatheter with 2 tip markers located 3 cm apart
- Guidewire compatible with the selected infusion microcatheter
- Three (3) continuous saline (or heparinized saline) flush set-ups with pressure bags: one (1) flush system for the femoral sheath, one (1) for the guiding catheter, and one (1) for the microcatheter
- Two (2) rotating hemostatic valves (RHV)
- Three-way stopcock
- One-way stopcock
- IV Pole (not required if DCB will be placed in sterile sleeve and placed in sterile field)

Note: A backup Codman Detachment Control Box is suggested for all procedures.

INSTRUCTIONS

MICROCOIL SIZE SELECTION

Microcoil selection is at the discretion of the physician. The appropriate microcoil size should be chosen based upon pre-embolization angiographic assessment of the diameter, height, and width of the aneurysm, as well as, the width of the aneurysm ostium. In most cases, the initial microcoil implanted should be a three-dimensional complex shape. To minimize the potential of microcoil migration away from the aneurysm, the diameter of the first microcoil selected should not be less than the width of the aneurysm neck ostium. The subsequent selected coils typically will be of decreasing size and the physician may continue to implant microcoils until he or she determines that the aneurysm has been successfully treated.

MICROCATHETER SELECTION

Proper selection of the appropriately sized microcatheter is required to avoid damage to the Microcoil System and to minimize potential complications. Microcatheter selection is also determined by the physician and is predicated by the location of the aneurysm, patient safety, and physician preference. To ensure proper placement and detachment of the Microcoil, the selected microcatheter must have two (2) radiopaque tip markers three (3) cm apart.

The XTRASOFT ORBIT GALAXY G2 System is compatible with microcatheters with inner lumen diameters ranging from 0.0165 to 0.019 in. (0.419 to 0.483 mm).

CONTINUOUS SALINE FLUSH SET-UP

To achieve optimal performance of the Microcoil System, it is important that a continuous infusion of an appropriate flush solution be maintained. Figure 2 illustrates the connections necessary for the Microcoil Delivery System including a typical continuous saline flush set-up with pressure bag for the catheter systems.

- Attach a rotating hemostatic valve (RHV) to the hub of the guiding catheter.
- Connect a three-way stopcock to the RHV sidearm, and attach a line to the stopcock for the continuous infusion of solution.
- Carefully select an appropriately sized infusion microcatheter based on the size of the selected Microcoil System.
- Insert the infusion microcatheter into the RHV connected to the guiding catheter hub.
- Attach a second RHV to the hub of the infusion microcatheter.
- Connect a one-way valve stopcock to the RHV's sidearm, and attach a flush line.
- Adjust the hydrostatic pressure to 300 mm of Hg and maintain an open flush throughout the procedure.

-

Fig. 2: Delivery System Set Up with Continuous Flush

Fig. 3: Holding the introducer tip and re-sheathing tool with the clear tab between your thumb and forefinger. The coil tip is also shown emerging from the introducer tip

Fig. 4: Aligning the Introducer Tip

Positioning the Microcoil

1. The proper alignment of the marker bands is illustrated in Figure 5. The infusion microcatheter has two distal tip marker bands located three (3) cm apart. The Microcoil System has a distal marker band three (3) mm from the detachment mechanism. To obtain proper positioning of the Microcoil System for detachment under fluoroscopic guidance, align the radiopaque marker on the DPU wire just past the proximal marker band on the tip of the microcatheter.
2. Once the desired microcoil placement has been achieved, gently tighten the RHV around the DPU wire to maintain its position.


Fig. 5: Proper Marker Band Alignment

DETACHING THE MICROCOIL

1. Re-verify the position of the microcoil within the aneurysm via fluoroscopy. The Microcoil System is now ready for detachment.
2. Press the power button on the Detachment Control Box (DCB), if not already powered on.
3. Attach the connecting cable to the connector end of the DPU wire, ensuring that it is fully seated. If not already connected, attach the other end of the connecting cable to the output connector of the DCB.
4. Verify that the Microcoil Delivery System is fully connected and no faults are indicated on the DCB. If a fault exists, reseal all connections between the DPU, the DCB, and the connecting cable. If a fault still persists, replace the connecting cable. If this does not correct the error, replace the DCB. If the Microcoil Delivery System still continues to have a fault, retrieve the microcoil as described in the following section, re-sheathing the Microcoil System, and replace with a new Microcoil System.
5. If the system is free of faults, depress the **Detach** button on the Blue EnPOWER DCB or the EnPOWER Control Cable. The **Detach Cycle** light next to the button will illuminate and an intermittent tone will sound for the duration of the detachment cycle. If the light and audible tone do not activate, replace the DCB.
6. After the light goes out and the tone stops, detachment of the microcoil from the DPU wire must be fluoroscopically verified by gently withdrawing the DPU wire back approximately one (1) mm. Observe if the microcoil has detached.
 - If the microcoil did not detach and no fault light is illuminated or flashing, repeat the steps above.
 - If a fault light is detected, the **System Ready** light is not illuminated, or there is no detachment after two detachment attempts, replace the DCB unit.
 - If detachment still does not occur, carefully withdraw the entire Microcoil System and redeploy a new Microcoil System.
7. After detaching the coil, remove the DPU wire from the microcatheter and discard.

Note: Do not disconnect the connecting cable from either the DCB or DPU wire during a detachment cycle.

CAUTION: Always perform fluoroscopic verification of detachment prior to withdrawing the DPU wire fully. Failure to do so may lead to an embolic complication.
8. Repeat the above sequence for all additional microcoils until the procedure is complete.

Note: If it is determined that the microcoil will not be detached, simply disconnect it from the connecting cable and remove it from the microcatheter as outlined in the following section.

RE-SHEATHING THE MICROCOIL SYSTEM

When necessary, the Microcoil System can be reloaded into the introducer sheath using the re-sheathing tool.

1. Loosen the RHV. Using the fluoroscopic guidance, retract the Microcoil from the aneurysm back into the Microcatheter.
2. Locate the introducer tip. Grasp the introducer tip with your left hand and the resheathing tool with your right hand. Pull with your right hand while holding on the resheathing tool towards the connector. This will start the resheathing of the coil and the DPU pusher wire.
3. When the resheathing tool reaches the end of the introducer sheath, replace the introducer tip back inside the RHV. Tighten the RHV.
4. With your right hand, grasp the connector and continue to pull the DPU wire out of the sheath until coil is completely inside the distal end of the introducer tip.
5. Loosen the RHV and remove the introducer.


Fig. 6: Pulling Back the DPU Wire Hub Connector

CAUTION: Pulling the exposed microcoil through the RHV grimmot may damage the coil.

CAUTION: Do not pull the DPU wire back too far since it may cause a softer section of the DPU wire to become exposed. An arm's length pull should re-sheath the coil.

6. Once the device is fully retracted, slide the re-sheathing tool over the DPU/ introducer sheath body until the system is relocked.
7. Loosen RHV and remove device. If needed, the Microcoil System is now ready to be re-inserted.

INSTRUCTIONS: BLUE EnPOWER DETACHMENT CONTROL BOX VERIFY PROPER FUNCTIONALITY OF DCB

Note: This DCB contains no user replaceable batteries.

Note: The EnPOWER DCB has an automatic shut off after ½ hour of non-use. If the DCB automatically shuts off, press the **Power** button to restart.

1. Press the **Power** button on the front panel of the DCB. All indicator lights will illuminate briefly then turn off while the unit performs a self test. With no cables connected, only one of the battery indicator lights should be illuminated.
 - If the green **Full Battery** light is illuminated, the DCB is functioning properly.
 - If the amber **Low Battery** light is illuminated, then the DCB will function properly, but less than 100 detach cycles remain on the unit.
 - If the red **Dead Battery** light is illuminated, the DCB will NOT function. Replace with new DCB.
 - If the red **System Fault** light, or no light, is illuminated, then the DCB is NOT functioning properly. The unit should not be used and should be returned to Codman & Shurtleff or an authorized representative for replacement.
2. Once verification of the DCB is complete and the DCB is functioning properly, mount the DCB onto an IV pole or equivalent mounting device.

Note: Using a flathead screwdriver or a thin coin, the clamp on the back of the DCB may be removed and the DCB inserted into a sterile sleeve so the device may be placed in the sterile field.

VERIFICATION OF CABLE CONNECTION

1. Open a sterile Codman connecting cable into the sterile field.
2. Pass the larger diameter connector out of the sterile field. The remainder of the connecting cable must remain in the sterile field.
3. Insert the larger cable connector into the DCB output connector until a click is felt.

VERIFICATION OF MICROCOIL SYSTEM

1. Select the appropriate Microcoil System based on the aneurysm's angiographic assessment.
2. Inspect the microcoil packaging for any breach of sterility.
3. Open the sterile peel-pouch containing the selected Microcoil System and remove the hoop containing the Microcoil System from the pouch.
4. With the Microcoil System remaining in the hoop, connect the smaller connector on the connecting cable to the hub connector of the Microcoil System's Device Positioning Unit (DPU) wire, and power on the DCB. The green **System Ready** light will illuminate when both the cable and microcoil are connected to the DCB and functioning properly.
 - If the **System Ready** light does not illuminate, check all cable connections.
 - If the **System Ready** light still does not illuminate, unplug the Microcoil System from the cable and select a new cable. Follow the steps outlined above and attach the new cable. Connect the existing Microcoil System to the new cable. If the **System Ready** light is now illuminated, the Microcoil Delivery System is functioning properly. Re-package the malfunctioning cable, and return it to Codman & Shurtleff or an authorized representative for replacement.

Note: The Codman connecting cable should not be re-sterilized or reused. Re-sterilization may damage the unit making it non-functional. Re-sterilized connecting cables will not be replaced by Codman & Shurtleff.
 - If the **System Ready** light still does not illuminate, select a replacement Microcoil System. Repeat the steps above to remove the microcoil from its sterile package. With the Microcoil System remaining in the

hoop, attach the new connecting cable. If the **System Ready** light is now illuminated, the Microcoil Delivery System is functioning properly. Re-package the malfunctioning Microcoil System, and return it to Codman & Shurtleff or an authorized representative for replacement.

- If the **System Ready** light still does not illuminate, the DCB is NOT functioning properly. The DCB should not be used and should be returned to Codman & Shurtleff or an authorized representative for replacement.

5. Disconnect the cable connector from the DPU hub connector. The Microcoil System and DCB are ready to be used.
6. Press the **Power** button again to turn off the DCB until it is ready to be used in the procedure.
Note: The EnPOWER DCB has an automatic shut off after ½ hour of non-use. If the DCB automatically shuts off, press the **Power** button to restart.
7. Refer back to the Microcoil Inspection and Placement Section of this IFU.

SAFETY INFORMATION:

The CODMAN Detachment Control Box (DCB000001-20 and DCB000005-00) comply with the requirements of:

- ANSI/AAMI ES60601-1:2005
- CAN/CSA-C22.2 No. 60601-1:2008
- IEC 60601-1 (2005) 3rd Edition
- EN 60601-1 (2006) 3rd Edition
- IEC 60601-1-6 (2010)
- IEC 62366 (2007)
- IEC 60601-1-2 (2007)

Detachment Control Box Specifications

Blue EnPOWER Detachment Control Box,
Catalog Number DCB000005-00

- Battery type: lithium
- Input voltage: 3.6 VDC
- Output voltage: 9.0 VDC, maximum
- Output current: 200 mA, maximum
- Operating temperature: 18 to 35°C (65 to 95°F)
- Operating humidity: 0 to 85% RH, non-condensing
- Transport & Storage Temperature: -18°C to 60°C
- Transport & Storage Humidity: 10% to 85% RH non-condensing
- Transport & Storage Pressure: 500 hPa to 1050 hPa
- Operating Pressure Range: 700 hPa to 1060 hPa
- Electrical safety classification: Class IIA Medical Product, internally powered, Type CF, ordinary equipment, intermittent operation
- Ingress of liquid classification: IPX2 rating

Cleaning Instructions for the Detachment Control Box

The DCB may be cleaned by wiping its surface with 70% isopropyl alcohol, a mild detergent, or disinfectant and a soft cloth. Use caution when cleaning near the output connector.

Repair of Detachment Control Box

There are no user serviceable parts in the DCB. If the DCB is not functioning properly it should not be used and should be returned to Codman & Shurtleff or an authorized representative for replacement.

Storage Instructions for the Detachment Control Box

Between uses, the DCB should be stored in a clean environment that is protected against extreme temperature and humidity.

Disposal of the Detachment Control Box

Do not discard this product as unsorted municipal waste. Contact your local authorities for disposal instructions or return it to Codman & Shurtleff for disposal.

MAGNETIC RESONANCE IMAGING (MRI) INFORMATION

Non-clinical testing demonstrated that the detachable coil is MR Conditional. A patient with this device can be scanned safely, immediately after implant, under the following conditions:

Static Magnetic Field

- static magnetic field of 3 Tesla or less
- maximum spatial gradient magnetic field of 720 Gauss/cm or less

MRI-RELATED HEATING

In non-clinical testing, the coil produced a temperature rise of +1.8 degrees centigrade during MRI performed for 15 minutes of scanning in a 3 Tesla MR system. Therefore, the MRI-related heating experiments for the detachable coil at 3 Tesla using a transmit/receive RF body coil at an MR system reported whole-body-averaged specific absorption rate (SAR) of 2.9-W/kg indicated that the greatest amount of heating that occurred in association with these specific conditions was equal to or less than +1.8 degrees centigrade. MR image quality can be compromised if the area of interest is in the same area or relatively close to the position of the detachable coil. Optimization of MR imaging parameters is recommended.

WARRANTY

Codman & Shurtleff, Inc. warrants that this medical device is free from defects in both materials and workmanship. **Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed. Suitability for use of this medical device for any particular surgical procedure should be determined by the user in conformance with the manufacturer's instructions for use. There are no warranties that extend beyond the description on the face hereof.**

® ORBIT GALAXY, XTRASOFT, and EnPOWER are registered trademarks of Codman & Shurtleff, Inc.

SYMBOLS

Blue EnPOWER DCB000005-00	
	Power on/off
	Detach button
	Detach cycle
	Full battery power
	Low battery power
	Dead battery
	System fault
	System ready
	Cardiac floating

**Table 1 – Guidance and Manufacturer’s Declaration – Emissions
All ME Equipment and ME Systems**

Guidance and Manufacturer’s Declaration - Emissions		
The EnPOWER Detachment Control Box is intended for use in the electromagnetic environment specified below. The customer or user of the EnPOWER DCB should ensure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment – Guidance
RF Emissions CISPR 11	Group 1	The EnPOWER DCB uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Group 2	The EnPOWER DCB must emit Electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.
RF Emissions CISPR 11	Class B	The EnPOWER DCB is suitable for use in all establishments, including domestic, and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonics IEC 61000-3-2	N/A	
Flicker IEC 61000-3-3	N/A	

**Table 2 – Guidance and Manufacturer’s Declaration – Immunity
All ME Equipment and ME Systems**

Guidance and Manufacturer’s Declaration – Immunity			
The EnPOWER DCB is intended for use in the electromagnetic environment specified below. The customer or user of the EnPOWER DCB should ensure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
ESD IEC 61000-4-2	±6 kV Contact ±8 kV Air	±6 kV Contact ±8 kV Air	Floors should be wood, concrete or ceramic tile. If floors are synthetic, the r/h should be at least 30%
EFT IEC 61000-4-4	±1 kV I/Os	±1 kV I/Os	
Power Frequency 50/60 Hz Magnetic Field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be that of a typical commercial or hospital environment.

**Table 3 – Guidance and Manufacturer’s Declaration – Immunity
ME Equipment and ME Systems that are NOT Life-supporting**

Guidance and Manufacturer’s Declaration – Immunity			
The EnPOWER Detachment Control Box is intended for use in the electromagnetic environment specified below. The customer or user of the EnPOWER DCB should ensure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	(V1)=3 Vrms (E1)=3 V/m	Portable and mobile communications equipment should be separated from the EnPOWER DCB by no less than the distances calculated/listed below: $D=(3.5/V1)\sqrt{P}$ 150kHz to 80MHz $D=(3.5/E1)\sqrt{P}$ 80 to 800 MHz $D=(7/E1)\sqrt{P}$ 800 MHz to 2.5 GHz where P is the max power in watts and D is the recommended separation distance in meters. Field strengths from fixed transmitters, as determined by an electromagnetic site survey, should be less than the compliance levels (V1 and E1). Interference may occur in the vicinity of equipment containing a transmitter.

Table 4 – Recommended Separation Distances between portable and mobile RF Communications equipment and the EnPOWER Detachment Control Box ME Equipment and ME Systems that are NOT Life-supporting

Recommended Separations Distances for the EnPOWER			
The EnPOWER DCB is intended for use in the electromagnetic environment in which radiated disturbances are controlled. The customer or user of the EnPOWER DCB can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF Communications Equipment and the EnPOWER DCB as recommended below, according to the maximum output power of the communications equipment.			
Max Output Power (Watts)	Separation (m) 150 kHz to 80 MHz $D=(3.5/V1)\sqrt{P}$	Separation (m) 80 to 800 MHz $D=(3.5/E1)\sqrt{P}$	Separation (m) 800 MHz to 2.5 GHz $D=(7/E1)\sqrt{P}$
0.01	0.11667	0.11667	0.23333
0.1	0.36894	0.36894	0.73785
1	1.1667	1.1667	2.3333
10	3.6894	3.6894	7.3785
100	11.667	11.667	23.333

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The following symbols pertain to the detachment control box only:

Les symboles suivants se rapportent uniquement à la commande de détachement :
Die folgenden Symbole gelten nur für den Schaltkasten zur Steuerung der Ablösung:
De volgende symbolen horen alleen bij de Detachment Control Box:
I simboli seguenti si riferiscono esclusivamente all'unità di controllo del distacco:
Los siguientes símbolos se refieren a la caja de control de separación únicamente:
Os seguintes símbolos dizem respeito apenas à caixa de controlo de desativação:
Følgende symboler vedrører alene DCB-boksen:
Följande symboler gäller endast för DCB:n:
Seuraavat symbolit koskevat vain irrotusohjainta:
Τα παρακάτω σύμβολα αφορούν μόνο το χειριστήριο αποσύνδεσης:
Az alábbi szimbólumok csak a leválobozra vonatkoznak:

Direct current

Courant continu
Gleichstrom
Gelijkstroom
Corrente continua
Corriente directa
Corrente directa
Direkte strøm
Likström
Tasavirta
Συνεχές ρεύμα
Egyenáram

Electrical and electronic equipment. Return waste to collection system or treatment and recycling facilities. Applicable in the EU. Follow decontamination instructions before return.

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Transport & Storage Temperature: -18°C to 60°C

Température de transport et de stockage : -18 °C à 60 °C
Temperatur bei Transport und Lagerung: -18 °C bis 60 °C
Temperatuur tijdens transport en opslag: -18 °C tot 60 °C
Temperatura di trasporto e immagazzinamento: Da -18 °C a 60 °C
Temperatura de transporte y almacenamiento: -18 °C a 60 °C
Temperatura de transporte e armazenamento: -18 °C a 60 °C
Temperatur ved transport og opbevaring: -18 °C til 60 °C
Temperatur vid transport och förvaring: -18 °C till 60 °C
Kuljetus ja säilytyslämpötila: -18...60 °C
Θερμοκρασία μεταφοράς και αποθήκευσης: -18 °C έως 60 °C
Szállítási és tárolási hőmérséklet: -18 °C-tól 60 °C-ig

Transport & Storage Humidity: 10% to 85% RH non-condensing
Humidité de transport et de stockage : 10 % à 85 % HR (sans condensation)
Luftfeuchte bei Transport und Lagerung: 10 bis 85 % rel. Luftfeuchte, nicht kondensierend
Luchtvochtigheid tijdens transport en opslag: 10 % tot 85 % RH niet-condenserend
Umidità di trasporto e immagazzinamento: dal 10 all'85 % di umidità relativa, senza condensa
Humedad de transporte y almacenamiento: 10 % a 85 % RH, sin condensación
Humidade de transporte e armazenamento: 10 % a 85 % HR sem condensação
Luftfugtighed ved transport og opbevaring: 10 % til 85 % relativ luftfugtighed, ikke-kondenserende
Luftfuktighet vid transport och förvaring: 10 % till 85 % RH, icke-kondenserande
Kuljetus- ja säilytyskosteus: 10–85 % suhteellinen kosteus, ei tiivistyvää
Υγρασία μεταφοράς και αποθήκευσης: 10 % έως 85 % σχετική υγρασία, χωρίς συμπύκνωση
Szállítási és tárolási páratartalom: 10–85 % relatív páratartalom, nem kondenzálódó

Transport & Storage Pressure: 500 hPa to 1050 hPa

Pression de transport et de stockage : 500 hPa à 1 050 hPa
Luftdruck bei Transport und Lagerung: 500 hPa bis 1050 hPa
Druk tijdens transport en opslag: 500 hPa tot 1050 hPa
Pressione di trasporto e immagazzinamento: da 500 hPa a 1050 hPa
Presión de transporte y almacenamiento: 500 hPa a 1050 hPa
Pressão de transporte e armazenamento: 500 hPa a 1050 hPa
Tryk ved transport og opbevaring: 500 hPa til 1050 hPa
Tryck vid transport och förvaring: 500 hPa till 1 050 hPa
Kuljetus- ja säilytyspaine: 500–1050 hPa
Πίεση μεταφοράς και αποθήκευσης: 500 hPa έως 1050 hPa
Szállítási és tárolási légnyomás: 500 hPa – 1050 hPa

Refer to Instructions Manual

Consulter le manuel d'instructions
Bedienungsanleitung beachten
Raadpleeg de gebruikershandleiding
Fare riferimento al manuale d'istruzioni.
Refírase al manual de instrucciones
Ver Manual de Instruções
Se instruktionshåndbogen
Se bruksanvisningen
Katso käyttöohjeet
Ανατρέξτε στο εγχειρίδιο οδηγιών
Lásd a Használati útmutatót

Do not discard this product as unsorted municipal waste. Contact your local authorities or return to Codman & Shurtleff for disposal.

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The following symbols pertain to the microcoil delivery system only:

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